



Cedara I-SoftView IHE Integration Statement		Date
		<i>October 1, 2004</i>
Vendor	Product Name	Version
<i>Cedara Software Corp.</i>	<i>Cedara I-ReadMammo™</i>	<i>5.2</i>
<i>This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options specified below:</i>		
Integration Profiles Implemented	Actors Implemented	Options Implemented
<ul style="list-style-type: none"> • Schedule Workflow (SWF) • Consistent Presentation of Images (CPI) • Access to Radiology Information (ARI) • Patient Information Reconciliation (PIR) • Key Image Note (KIN) • Evidence Documents (ED) 	Image Display	None
Links to Standards Conformance Statements of the Implementation		
HL7		
DICOM		
Link to further information on this product:	http://www.cedara.com/products/diagnostic_and_review_workstations.htm	