



<b>Cedara I-SoftView IHE Integration Statement</b>		<b>Date</b>
		<i>October 15, 2003</i>
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>
<i>Cedara Software Corp.</i>	<i>Cedara I-SoftView (Report)</i>	<i>5.2</i>
<i>This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options specified below:</i>		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
<ul style="list-style-type: none"> <li>• Schedule Workflow (SWF)</li> <li>• Consistent Presentation of Images (CPI)</li> <li>• Access to Radiology Information (ARI)</li> <li>• Patient Information Reconciliation (PIR)</li> <li>• Key Image Note (KIN)</li> <li>• Evidence Documents (ED)</li> </ul>	Image Display	None
<b><i>Links to Standards Conformance Statements of the Implementation</i></b>		
<b>HL7</b>		
<b>DICOM</b>		
<b><i>Link to further information on this product:</i></b>	<a href="http://www.cedara.com/products/diagnostic_and_review_workstations.htm">http://www.cedara.com/products/diagnostic_and_review_workstations.htm</a>	